

**Overview**
**Method Name**

ELISA

**NY State Available**

Yes

**Specimen**
**Specimen Type**

Serum

**Specimen Required**
**Collection Container/Tube:** 5 mL Red

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 1 mL

**Acceptable:** SST

**Collection Instructions:** Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

**Specimen Minimum Volume**

0.3 mL (volume does NOT allow for repeat testing)

**Reject Due To**

Gross hemolysis:	Reject
Thawing:	Warm OK; Cold OK
Gross lipemia:	Reject
Gross icterus	Reject
Other:	Anything other than serum; bacterial contamination

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Frozen	60 days	
	Ambient	7 days	

**Clinical and Interpretive**

## Clinical Information

The anti-PM/Scl-100 antibody is associated with younger age, calcinosis and has lower rates of gastrointestinal symptoms, ILD and pulmonary hypertension. There is also evidence of a possibly better survival compared to the presence of either anti-PM/Scl-75 or anti-Scl-70 antibodies.

## Reference Values

Reference Range: <20

## Interpretation:

Negative: <20 units

Weak Positive: 20-39 units

Moderate Positive: 40-80 units

Strong Positive: >80 units

## Performance

### PDF Report

No

### Performing Laboratory Location

Esoterix Endocrinology

## Fees and Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

### Test Classification

This test was developed and its performance characteristics determined by LabCorp. It has not been cleared or approved by the Food and Drug Administration.

### CPT Code Information

83520

### LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FAPMA	Anti-PM/Scl-100 Ab	Not Provided

---

Result ID	Test Result Name	Result LOINC Value
FAPMA	Anti-PM/ScI-100 Ab	Not Provided