

**Overview**
**Method Name**

Blocking Radioimmunoassay (RIA) after acid-alcohol extraction

**NY State Available**

Yes

**Specimen**
**Specimen Type**

Serum

**Specimen Required**

Draw blood in a plain, red-top tube(s). Separate within 1 hour of collection, freeze immediately. Send 0.5 mL serum frozen.

Note: 1. Serum gel tube is okay, but must pour off into a plastic screw cap vial and freeze.

2. Minimum volume does not permit for repeat analysis

**Specimen Minimum Volume**

0.1 mL

NOTE: Minimum volume does not allow for repeat analysis.

**Reject Due To**

Hemolysis	NA
Lipemia	NA
Icteric	NA
Other	NA

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Frozen	200 days	

**Clinical and Interpretive**
**Reference Values**

Age	Range (ng/mL)	Mean
Prepubertal	258 – 882	570

Pubertal	273 Åç&#130;-â&#128;&#156; 892	583
Adults	333 Åç&#130;-â&#128;&#156; 967	650

## Performance

### PDF Report

No

### Day(s) Performed

Tuesday

### Report Available

3 to 14 days

### Performing Laboratory Location

Esoterix Endocrinology

## Fees and Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

### Test Classification

This test was developed and its performance characteristics determined by LabCorp. It has not been cleared or approved by the Food and Drug Administration.

### CPT Code Information

83519

### LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FIGF2	Insulin-Like Growth Factor II	2485-1

Result ID	Test Result Name	Result LOINC Value
Z2324	Insulin-Like Growth Factor II	2485-1