

**Overview****Method Name**

Spectrophotometry

**NY State Available**

Yes

**Specimen****Specimen Type**

Urine

**Specimen Required**

Container/Tube: Plastic urine container

Specimen Volume: 2 mL

Collection Instructions:

1. Collect random urine without preservative.
2. Ship 2 mL urine in a sterile screw capped plastic container ship ambient.

**Specimen Minimum Volume**

1 mL

**Reject Due To**

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	Preservatives

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Urine	Ambient (preferred)	5 days	
	Frozen	21 days	
	Refrigerated	5 days	

**Clinical and Interpretive****Reference Values**

Adult: &lt; or = 4 U/L

**Performance****PDF Report**

No

**Day(s) Performed**

Tuesday - Friday, Sunday

**Report Available**

1 to 7 days

**Performing Laboratory Location**

Quest Diagnostics Nichols Institute

**Fees and Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

**CPT Code Information**

83690

**LOINC® Information**

Test ID	Test Order Name	Order LOINC Value
FLIPR	Lipase, Urine	10888-6

Result ID	Test Result Name	Result LOINC Value
Z1559	Lipase, Urine	10888-6