

Overview

Method Name

Enzyme Immunoassay

NY State Available

Yes

Specimen

Specimen Type

Serum

Specimen Required

Specimen Type: Serum

Container/Tube: Red or SST

Specimen Volume: 0.5 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume

0.5 mL

Reject Due To

|         |        |
|---------|--------|
| Lipemia | Reject |
|---------|--------|

Specimen Stability Information

| Specimen Type | Temperature              | Time     | Special Container |
|---------------|--------------------------|----------|-------------------|
| Serum         | Refrigerated (preferred) | 28 days  |                   |
|               | Ambient                  | 28 days  |                   |
|               | Frozen                   | 365 days |                   |

Clinical & Interpretive

Reference Values

Negative

Interpretation

The Phadia Phadiatop test is an allergy screening test with excellent sensitivity and specificity for inhalant allergy. It uses an ImmunoCAP with a balanced mixture of representative allergens, including grasses, trees, weeds, cat, dog, mites and molds. A positive result indicates that the patient is allergic to one or more of these allergens; a negative indicates the patient is not allergic to inhalant allergens. Note that the test does not assess a patient's sensitivity to food, drug, chemical or certain unusual or rare allergens.

Performance

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

3 to 6 days

Performing Laboratory Location

Eurofins Viracor

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

86005

LOINC® Information

| Test ID | Test Order Name            | Order LOINC® Value |
|---------|----------------------------|--------------------|
| FPHAS   | Phadiatop (Allergy Screen) | 57036-6            |

| Result ID | Test Result Name           | Result LOINC® Value |
|-----------|----------------------------|---------------------|
| FPHAS     | Phadiatop (Allergy Screen) | 57036-6             |