



Test Definition: FPHAS

Phadiatop (Allergy Screen)

Overview

Method Name

Enzyme Immunoassay

NY State Available

Yes

Specimen

Specimen Type

Serum

Specimen Required

Specimen Type: Serum

Container/Tube: Red or SST

Specimen Volume: 0.5 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume

0.5 mL

Reject Due To

Lipemia	Reject
---------	--------

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	365 days	

Clinical & Interpretive

Reference Values

Negative

Interpretation

The Phadia Phadiatop test is an allergy screening test with excellent sensitivity and specificity for inhalant allergy. It uses an ImmunoCAP with a balanced mixture of representative allergens, including grasses, trees, weeds, cat, dog, mites and molds. A positive result indicates that the patient is allergic to one or more of these allergens; a negative indicates the patient is not allergic to inhalant allergens. Note that the test does not assess a patient's sensitivity to food, drug, chemical or certain unusual or rare allergens.

Performance**PDF Report**

No

Day(s) Performed

Monday through Friday

Report Available

3 to 6 days

Performing Laboratory Location

Eurofins Viracor

Fees & Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

86005

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FPHAS	Phadiatop (Allergy Screen)	57036-6

Result ID	Test Result Name	Result LOINC® Value
FPHAS	Phadiatop (Allergy Screen)	57036-6