

Overview

Method Name
Multiplex array electrochemiluminescence

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Specimen Type: Serum
Collection Container/Tube: Red or SST
Submission Container/Tube: Plastic vial
Specimen Volume: 1 mL
Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1 mL of serum frozen in a plastic vial.

Specimen Minimum Volume
1 mL

Reject Due To

| | |
|----------|----------------------|
| Thawing: | Warm reject; Cold OK |
|----------|----------------------|

Specimen Stability Information

| Specimen Type | Temperature | Time | Special Container |
|---------------|-------------|----------|-------------------|
| Serum | Frozen | 365 days | |

Clinical & Interpretive

Reference Values
<8.6 pg/mL

Performance

PDF Report

No

Day(s) Performed

Monday, Wednesday, Friday

Report Available

5 to 11 days

Performing Laboratory Location

Eurofins Viracor

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by Viracor Eurofins. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

83520

LOINC® Information

| Test ID | Test Order Name | Order LOINC® Value |
|---------|-----------------|--------------------|
| FIFNY | IFN-γ, Serum | 27415-9 |

| Result ID | Test Result Name | Result LOINC® Value |
|-----------|------------------|---------------------|
| FIFNY | IFN-γ, Serum | 27415-9 |