



# Test Definition: FZOLP

Zolpidem (Ambien), serum or plasma

## Overview

### Method Name

Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS)

### NY State Available

Yes

## Specimen

### Specimen Type

Varies

### Specimen Required

Submit only 1 of the following specimens:

#### Plasma

Draw blood in a green-top (sodium heparin) tube(s), plasma gel tube is not acceptable. Spin down and send 2 mL sodium heparin plasma refrigerated in a plastic vial.

#### Serum

Draw blood in a plain red-top tube(s), serum gel tube is not acceptable. Spin down and send 2 mL of serum refrigerated in a plastic vial.

### Specimen Minimum Volume

0.25 mL

### Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	NA

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	7 days	
	Ambient	72 hours	
	Frozen	180 days	

## Clinical & Interpretive

### Reference Values

Units: ng/mL

Expected hypnotic zolpidem concentrations in patients taking recommended daily dosages: up to 250 ng/mL.

Toxic range has not been established.

## Performance

### PDF Report

No

### Day(s) Performed

Monday through Sunday

### Report Available

3 to 7 days

### Performing Laboratory Location

Medtox Laboratories, Inc.

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### CPT Code Information

80368

### LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FZOLP	Zolpidem (Ambien), Serum	6953-4

Result ID	Test Result Name	Result LOINC® Value
Z3322	Zolpidem, SE	6953-4