

Overview

Method Name

Enzymatic/Spectrophotometric

NY State Available

Yes

Specimen

Specimen Type

Fecal

Specimen Required

Specimen Type: Random stool

Container/Tube: Sterile leak proof stool container

Specimen Volume: 1 g

Collection Instructions: Collect 1 gm random stool in sterile leak proof container, ship refrigerate.

Note:

Dietary restrictions: Patients receiving pancreatic enzymes should discontinue taking the enzymes at least 5 days before the collection of the stool sample.

Specimen Minimum Volume

0.5 gram

Reject Due To

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Fecal	Refrigerated (preferred)	14 days	
	Frozen	30 days	

Clinical & Interpretive

Reference Values

2.3 – 51.4 U/g

Performance

PDF Report

No

Day(s) Performed

Wednesday

Report Available

5 to 18 days

Performing Laboratory Location

Quest Diagnostics Nichols Institute

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its analytical performance characteristics have been determined by Quest Diagnostics. It has not been cleared or approved by FDA. This assay has been validated pursuant to the CLIA regulations and is used for clinical purposes.

CPT Code Information

84311

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FCHYS	Chymotrypsin, Stool	25375-7

Result ID	Test Result Name	Result LOINC® Value
FCHYS	Chymotrypsin, Stool	25375-7