



Test Definition: FMACR

Macroamylase

Overview

Method Name

Polyethylene Glycol Precipitation, Spectrophotometry

NY State Available

Yes

Specimen

Specimen Type

Serum

Specimen Required

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 2 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume

1 mL

Reject Due To

Hemolysis	Mild OK; Gross reject
Lipemia	NA
Icterus	NA
Other	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	7 days	
	Frozen	30 days	

Clinical & Interpretive

Reference Values

Not Detected

Performance**PDF Report**

No

Day(s) Performed

Tuesday, Thursday

Report Available

2 to 9 days

Performing Laboratory Location

Quest Diagnostics Nichols Institute

Fees & Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

82150

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FMACR	Macroamylase	15358-5

Result ID	Test Result Name	Result LOINC® Value
FMACR	Macroamylase	15358-5