

Overview

Method Name

Gas Chromatography/Mass Spectrometry (GC/MS)

NY State Available

Yes

Specimen

Specimen Type

Urine

Specimen Required

Container/Tube: Plastic, preservative-free urine container

Specimen Volume: 20 mL

Collection Instructions:

- 1. Collect 20 mL random urine without preservative.
- 2. Transfer into a preservative-free plastic urine container.
- 3. Ship refrigerate.

Specimen Minimum Volume

10 mL

Reject Due To

All specimens will be evaluated by the processing and performing laboratories for test suitability.

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	10 days	
	Ambient	72 hours	
	Frozen	396 days	

Clinical & Interpretive

Reference Values

Results in ng/mL are reported to the limit of quantitation for the analysis.

Performance

PDF Report

No

Day(s) Performed

Varies

Report Available

16 to 20 days

Performing Laboratory Location

Medtox Laboratories, Inc.

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

80375

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FMETH	Methaqualone, Confirmation	16244-6

Result ID	Test Result Name	Result LOINC® Value
Z4116	Methaqualone, Urine Quant	16244-6