

Test Definition: FMARC

Bupivacaine (Marcaine)

Overview

Method Name

Liquid Chromatography/ Tandem Mass Spectrometry (LC/MS/MS)

NY State Available

Yes

Specimen

Specimen Type

Serum Red

Specimen Required

Draw blood in a plain red-top tube(s), serum gel tube is not acceptable. Spin down and send 1 mL of serum frozen in a plastic vial.

Specimen Minimum Volume

0.5 mL

Reject Due To

| Hemolysis | NA |
|-----------|----|
| Lipemia | NA |
| Icterus | NA |
| Other | NA |

Specimen Stability Information

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------------|----------|-------------------|
| Serum Red | Frozen (preferred) | 180 days | |
| | Ambient | 72 hours | |
| | Refrigerated | 72 hours | |

Clinical & Interpretive

Reference Values

Reference Range: 0.1 - 4.0 ug/mL



Test Definition: FMARC

Bupivacaine (Marcaine)

Performance

PDF Report

No

Day(s) Performed

Monday through Sunday

Report Available

5 to 11 days

Performing Laboratory Location

Medtox Laboratories, Inc.

Fees & Codes

Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

Test Classification

This test was developed and its performance characteristics determined by LabCorp. It has not been cleared or approved by the Food and Drug Administration.

CPT Code Information

80299

LOINC® Information

| Test ID | Test Order Name | Order LOINC® Value |
|---------|------------------------|--------------------|
| FMARC | Bupivacaine (Marcaine) | 3412-4 |

| Result ID | Test Result Name | Result LOINC® Value |
|-----------|------------------|---------------------|
| Z5052 | Bupivacaine | 3412-4 |