



# Test Definition: FDGTX

Digitoxin, Serum

## Overview

### Method Name

Immunoassay (IA)

### NY State Available

Yes

## Specimen

### Specimen Type

Serum Red

### Specimen Required

**\*\*\*Must submit one specimen per order. Specimens cannot be shared between multiple orders.\*\*\***

**Specimen Type:** Serum

**Collection Container/Tube:** Red top

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 1 mL

#### Collection Instructions:

1. Draw blood in a plain, red-top tube(s). **Serum gel tube is not acceptable.**
2. Centrifuge and send 1 mL of serum refrigerated in a plastic, preservative-free vial.

### Specimen Minimum Volume

0.4 mL

### Reject Due To

Other	Room Temperature; Serum gel tube (SST), Plasma gel tube (PST)
-------	---

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	7 days	
	Frozen	90 days	

## Clinical & Interpretive

### Reference Values

Reporting limit determined each analysis.

---

None Detected ng/mL

Therapeutic Range: 10 - 30 ng/mL

## Performance

### PDF Report

Referral

### Day(s) Performed

Monday - Sunday

### Report Available

9 to 11 days

### Performing Laboratory Location

NMS Labs

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### CPT Code Information

80299

### LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FDGTX	Digitoxin	3559-2

Result ID	Test Result Name	Result LOINC® Value
FDGTX	Digitoxin	3559-2