

Overview

Method Name

Immunoassay (IA)

NY State Available

Yes

Specimen

Specimen Type

Serum Red

Specimen Required

Must submit one specimen per order. Specimens cannot be shared between multiple orders.

Specimen Type: Serum

Collection Container/Tube: Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL

Collection Instructions:

1. Draw blood in a plain, red-top tube(s). Serum gel tube is not acceptable.

2. Centrifuge and send 1 mL of serum refrigerated in a plastic, preservative-free vial.

Specimen Minimum Volume

0.4 mL

Reject Due To

Other	Room Temperature; Serum gel tube (SST), Plasma gel tube (PST)
-------	---

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	7 days	
	Frozen	90 days	

Clinical & Interpretive

Reference Values

Reporting limit determined each analysis.

None Detected ng/mL
Therapeutic Range: 10 - 30 ng/mL

Performance

PDF Report
Referral

Day(s) Performed
Monday - Sunday

Report Available
9 to 11 days

Performing Laboratory Location
NMS Labs

Fees & Codes

- Fees
- Authorized users can sign in to [Test Prices](#) for detailed fee information.
 - Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
 - Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information
80299

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FDGTX	Digitoxin	3559-2

Result ID	Test Result Name	Result LOINC® Value
FDGTX	Digitoxin	3559-2