

Overview

Method Name
Reversed-phase UPLC-MS/MS

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Specimen Type: Serum
Container/Tube: Red
Specimen Volume: 1 mL
Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube is not acceptable. Allow to clot for 30 to 60 minutes, spin down and send 1 mL of serum Frozen in a sterile, screw top tube.

Specimen Minimum Volume
0.5 mL

Reject Due To

| | |
|------------------|----------------------|
| Gross hemolysis: | Reject |
| Thawing: | Warm Reject; Cold OK |
| Gross lipemia: | Reject |
| Other: | SST |

Specimen Stability Information

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------------|---------|-------------------|
| Serum | Frozen (preferred) | 30 days | |
| | Ambient | 14 days | |
| | Refrigerated | 14 days | |

Clinical & Interpretive

Reference Values

The analytical measuring range is 0.1-10 mcg/mL. The therapeutic range has not yet been established.

Performance

PDF Report

No

Performing Laboratory Location

Eurofins Viracor

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

80299

LOINC® Information

| Test ID | Test Order Name | Order LOINC® Value |
|---------|-----------------|--------------------|
| FIVCZ | Isavuconazole | Not Provided |

| Result ID | Test Result Name | Result LOINC® Value |
|-----------|------------------|---------------------|
| FIVCZ | Isavuconazole | Not Provided |