

# **Test Definition: FVENE**

Venison IgE

## **Overview**

#### **Method Name**

Enzyme Immunoassay (EIA)

#### **NY State Available**

Yes

## **Specimen**

## **Specimen Type**

Serum

## **Specimen Required**

Specimen Type: Serum Container/Tube: Red Top Specimen Volume: 0.5 mL

Collection Instructions: Draw blood in a plain re-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL

of serum refrigerated in a plastic vial.

## **Specimen Minimum Volume**

0.5 mL

### **Reject Due To**

Lipemia	Reject

## **Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	365 days	

## **Clinical & Interpretive**

#### **Reference Values**

<0.35 kU/L

#### Interpretation



# **Test Definition: FVENE**

Venison IgE

<u>Class</u>	IgE (kU/L)	<u>Comment</u>
0	<0.35	<b>Below Detection</b>
1	0.35-0.69	Low Positive
2	0.70-3.49	Moderate Positive
3	3.50-17.49	Positive
4	17.50-49.99	Strong Positive
5	>49.99	Very Strong Positive

### **Performance**

### **PDF Report**

No

## **Performing Laboratory Location**

**Eurofins Viracor** 

#### **Fees & Codes**

#### **Fees**

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

## **Test Classification**

This test was developed and its performance characteristics determined by Eurofins Viracor . It has not been cleared or approved by the U.S. Food and Drug Administration.

### **CPT Code Information**

86003

#### **LOINC®** Information

Test ID	Test Order Name	Order LOINC® Value
FVENE	Venison IgE	7759-4

Result ID	Test Result Name	Result LOINC® Value
Z5660	Venison IgE	7759-4
Z5661	Class	102840-6