

Overview

Method Name  
Immunofixation Electrophoresis

NY State Available  
Yes

Specimen

Specimen Type  
CSF

Specimen Required  
Specimen Type: CSF  
Container/Tube: sterile screw cap container  
Specimen Volume: 6 mL  
Collection Instructions: Collect 6 mL of spinal fluid (CSF) in a sterile screw cap container. Ship frozen.

Specimen Minimum Volume  
4.5 mL

Reject Due To

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
CSF	Frozen (preferred)	30 days	
	Ambient	6 days	
	Refrigerated	10 days	

Clinical & Interpretive

Reference Values  
No abnormal bands are present on immunofixation.

Performance

PDF Report

Referral

Day(s) Performed

Sunday through Friday

Report Available

5 to 9 days

Performing Laboratory Location

Quest Diagnostics Nichols Institute

Fees & Codes

- Fees
- Authorized users can sign in to [Test Prices](#) for detailed fee information.
  - Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
  - Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

86335

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FIXCF	Immunofixation, CSF	13174-8

Result ID	Test Result Name	Result LOINC® Value
FIXCF	Immunofixation, CSF	13174-8