

Overview

Method Name
Immunofixation Electrophoresis

NY State Available
Yes

Specimen

Specimen Type
CSF

Specimen Required
Specimen Type: CSF
Container/Tube: sterile screw cap container
Specimen Volume: 6 mL
Collection Instructions: Collect 6 mL of spinal fluid (CSF) in a sterile screw cap container. Ship frozen.

Specimen Minimum Volume
4.5 mL

Reject Due To

Specimen Stability Information

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------------|---------|-------------------|
| CSF | Frozen (preferred) | 30 days | |
| | Ambient | 6 days | |
| | Refrigerated | 10 days | |

Clinical & Interpretive

Reference Values
No abnormal bands are present on immunofixation.

Performance

PDF Report

Referral

Day(s) Performed

Sunday through Friday

Report Available

5 to 9 days

Performing Laboratory Location

Quest Diagnostics Nichols Institute

Fees & Codes

- Fees
- Authorized users can sign in to [Test Prices](#) for detailed fee information.
 - Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
 - Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

86335

LOINC® Information

| Test ID | Test Order Name | Order LOINC® Value |
|---------|---------------------|--------------------|
| FIXCF | Immunofixation, CSF | 13174-8 |

| Result ID | Test Result Name | Result LOINC® Value |
|-----------|---------------------|---------------------|
| FIXCF | Immunofixation, CSF | 13174-8 |