

Acetylcholine Receptor (Muscle AChR) Binding
Antibody, Serum

Overview

Useful For

Supporting the diagnosis of autoimmune myasthenia gravis (MG) in adults and children

Distinguishing autoimmune from congenital MG in adults and children or other acquired forms of neuromuscular junction transmission disorders

An adjunct to the test for P/Q-type calcium channel binding antibodies as a diagnostic aid for Lambert-Eaton myasthenic syndrome

Testing Algorithm

This is the primary diagnostic test for myasthenia gravis.

Method Name

Radioimmunoassay (RIA)

NY State Available

Yes

Specimen

Specimen Type

Serum

Ordering Guidance

Standalone testing (this test) is recommended in certain situations.

This test should not be requested in patients who have recently received radioisotopes, therapeutically or diagnostically, because of potential assay interference. The specific waiting period before specimen collection will depend on the isotope administered, the dose given, and the clearance rate in the individual patient. Specimens will be screened for radioactivity prior to analysis. Radioactive specimens received in the laboratory will be held for 1 week and assayed if sufficiently decayed or canceled if radioactivity remains.

Specimen Required

Patient Preparation: For optimal antibody detection, specimen collection is recommended prior to initiation of immunosuppressant medication.

Supplies: Sarstedt Aliquot Tube 5 mL (T914)

Collection Container/Tube:

Preferred: Red top



Acetylcholine Receptor (Muscle AChR) Binding Antibody, Serum

Acceptable: Serum gel

Submission Container/Tube: Plastic vial

Specimen Volume: 1.5 mL

Forms

If not ordering electronically, complete, print, and send a <u>Neurology Specialty Testing Client Test Request</u> (T732) with the specimen.

Specimen Minimum Volume

1 mL

Reject Due To

Gross	Reject
hemolysis	
Gross lipemia	Reject
Gross icterus	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	72 hours	
	Frozen	28 days	

Clinical & Interpretive

Clinical Information

Fatigable weakness due to impaired postsynaptic transmission at the neuromuscular junction is characteristic of myasthenia gravis (MG). A clinical diagnosis should be supported by electrodiagnostic testing (ie, clinical-electrodiagnosis [EDX]). Positive autoimmune serology increases certainty of MG diagnosis but needs to be interpreted in the proper clinical-EDX context with response to anticholinesterase medications supporting the diagnosis. Most cases are autoimmune and are caused by IgG autoantibodies binding to critical postsynaptic membrane molecules (nicotinic muscle acetylcholine receptor [AChR] or its interacting proteins, such as muscle-specific kinase). Serologically, the detection of AChR binding antibody provides the best diagnostic sensitivity. However, the presence of both AChR binding and modulating activity improves diagnostic accuracy. Autoantibody detection frequency is lowest in patients with weakness confined to extraocular muscles (approximately 70% are positive for AChR binding antibodies) and highest in patients with generalized weakness due to MG (approximately 90% are positive for AChR binding antibodies). In adults with MG and AChR antibodies, approximately 20% will have thymoma and, very rarely (<1%), extrathymic cancers. Computed tomography imaging of the chest is considered the standard of care to evaluate for thymoma.

These results should only be interpreted in the appropriate clinical and electrophysiological context and are not diagnostic in isolation.



Acetylcholine Receptor (Muscle AChR) Binding Antibody, Serum

Note: Single antibody tests may be requested in the follow-up of patients with positive results previously documented in this laboratory.

Reference Values

< or =0.02 nmol/L

Interpretation

Positive results (>0.02 nmol/L) are indicative of autoimmune myasthenia gravis (MG). These results should be interpreted in the appropriate clinical and electrophysiological context.

With a diagnosis of MG, a paraneoplastic basis should be considered with thymoma being the most frequently associated tumor with MG.

The clinical sensitivity of this assay is approximately 90% in nonimmunosuppressed patients with generalized MG. The frequency of antibody detection is lower in MG patients with weakness clinically restricted to ocular muscles (71%), and antibody titers are generally low in ocular MG (eg, 0.03-1.0 nmol/L).

Negative results do not exclude the diagnosis of MG. If clinical suspicion remains and symptoms persist or worsen consider retesting. Results may be negative in the first 12 months after symptoms of MG appear or during immunosuppressant therapy. **Note:** In follow up of seronegative patients with adult-acquired generalized MG, 17.4% seroconvert to positive at 12 months (ie, seronegativity rate at 12 months is 8.4%). A subset of MG patients that are persistently negative for acetylcholine receptor binding antibodies will have muscle-specific kinase (MuSK) antibodies, and therefore, it is recommended to test for MuSK antibodies in seronegative patients with high clinical suspicion of MG.

In general, there is not a close correlation between antibody titer and severity of weakness, but in individual patients, clinical improvement may be accompanied by a decrease in titer.

Cautions

The presence of elevated immunoglobulins due to therapeutic intervention or other disorders (ie, hypergammaglobulinemia) may lead to false-positive results.

Positive results may be found in some patients with Lambert-Eaton syndrome, paraneoplastic central nervous system, and peripheral nervous system autoimmune disorders and in healthy individuals.

The presence of alpha-bungarotoxin antibodies may interfere with this assay.

Specimens ideally should be collected prior to initiation of immunosuppressive therapies as these may reduce the sensitivity of this test.

Clinical Reference

- 1. Lennon VA. Serological profile of myasthenia gravis and distinction from the Lambert-Eaton myasthenic syndrome. Neurology. 1997;48(Suppl 5):S23-S27. doi:10.1212/WNL.48.Suppl_5.23S
- 2. Lachance DH, Lennon VA. Paraneoplastic neurological autoimmunity. In: Kalman B, Brannagan III T, eds. Neuroimmunology in Clinical Practice. Blackwell Publishing Ltd; 2008:210-217



Acetylcholine Receptor (Muscle AChR) Binding
Antibody, Serum

- 3. Gilhus NE. Myasthenia gravis. N Engl J Med. 2016;375(26):2570-2581. doi:10.1056/NEJMra1602678
- 4. Nicolle MW. Myasthenia gravis and Lambert-Eaton myasthenic syndrome. Continuum (Minneap Minn). 2016;22(6, Muscle and Neuromuscular Junction Disorders):1978-2005. doi:10.1212/CON.000000000000415
- 5. Shelly S, Paul P, Bi H, et al. Improving accuracy of myasthenia gravis autoantibody testing by reflex algorithm. Neurology. 2020;95(22):e3002-e3011. doi:10.1212/WNL.00000000010910

Performance

Method Description

Fetal and adult detergent-solubilized, acetylcholine receptors (extracted from cultures of rhabdomyosarcoma [RD] cells) labeled with (125)I-alpha-bungarotoxin are incubated with patient sample. Anti-human IgG is then added to form an immunoprecipitate. After washing the precipitated immune complexes, the amount of (125)I-labeled receptor in the immunoprecipitate is measured using a gamma-counter. The amount of gamma emission in the precipitate is proportional to the amount of AChR-IgG in the sample. Results are reported in units of precipitated antigen (nMoI) per L of patient sample. (Griesmann GE, Kryzer TJ, Lennon VA. Autoantibody profiles of myasthenia gravis and Lambert-Eaton myasthenic syndrome. In: Rose NR, Hamilton RG, eds. Manual of Clinical and Laboratory Immunology. 6th ed. ASM Press; 2002:1005-1012; Waters P, Pettingill P, Lang B. Detection methods for neural autoantibodies. Handb Clin Neurol. 2016;133:147-163. doi:10.1016/B978-0-444-63432-0.00009-8)

PDF Report

No

Day(s) Performed

Monday through Sunday

Report Available

3 to 6 days

Specimen Retention Time

28 days

Performing Laboratory Location

Mayo Clinic Laboratories - Rochester Main Campus

Fees & Codes

Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.



Acetylcholine Receptor (Muscle AChR) Binding Antibody, Serum

Test Classification

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. It has not been cleared or approved by the US Food and Drug Administration.

CPT Code Information

86041

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
ARBI	ACh Receptor (Muscle) Binding Ab	97558-1

	Result ID	Test Result Name	Result LOINC® Value
ſ	8338	ACh Receptor (Muscle) Binding Ab	97558-1