

# **Test Definition: FPROP**

Propofol, Serum/Plasma

Overview

## Method Name

Gas Chromatography (GC)

### NY State Available

Yes

Specimen

Specimen Type Varies

### **Specimen Required**

\*\*\*Must submit one specimen per order. Specimens cannot be shared between multiple orders.\*\*\* Submit only 1 of the following specimens:

### Serum

Draw blood in a plain, red-top tube(s). **Serum gel tube is not acceptable.** Centrifuge and send 2 mL of serum refrigerate in plastic vial.

Note: 1. Indicate serum on request form.

2. Label specimen appropriately (serum).

### Plasma

Draw blood in a lavender-top (EDTA) tube(s). Plasma gel tube is not acceptable. Centrifuge and send 2 mL of EDTA plasma refrigerate in plastic vial.

Note: 1. Indicate plasma on request form.

2. Label specimen appropriately (plasma).

### Specimen Minimum Volume

1.0 mL

### **Reject Due To**

Hemolysis	NA
Lipemia	NA
Icteric	NA
Other	Polymer gel separation tube (SST or PST)

## Specimen Stability Information



# **Test Definition: FPROP**

Propofol, Serum/Plasma

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated		

### Clinical & Interpretive

### **Reference Values**

Reporting limit determined each analysis

Synonym(s): Diprivan

Patients required a mean blood propofol concentration of 4.05 +/- 1.01 mcg/mL for major surgery and 2.97 +/- 1.07 mcg/mL for non-major surgery. Blood propofol concentrations at which 50% of patients were awake and oriented after surgery were 1.07 and 0.95 mcg/mL respectively.

Psychomotor performance returned to baseline at blood propofol concentrations of 0.38 - 0.43 mcg/mL.

### Performance

PDF Report No

Day(s) Performed Monday through Sunday

Report Available 8 to 12 days

Specimen Retention Time 2 weeks

Performing Laboratory Location NMS Labs

### Fees & Codes

#### Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.



### **Test Classification**

This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the U.S. Food and Drug Administration.

# **CPT Code Information**

80299

### LOINC<sup>®</sup> Information

Test ID	Test Order Name	Order LOINC <sup>®</sup> Value
FPROP	Propofol	34673-4
Result ID	Test Result Name	Result LOINC <sup>®</sup> Value
Result IDZ1448	Test Result NamePropofol	Result LOINC <sup>®</sup> Value 34673-4