

## Overview

### Method Name

Enzymatic, Colorimetric Assay

### NY State Available

Yes

## Specimen

### Specimen Type

Varies

### Specimen Required

Submit only 1 of the following specimens:

#### Serum

Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 1 mL of serum ambient in a plastic vial.

#### Plasma

Draw blood in a lavender-top (EDTA) tube(s). Spin down and send 1 mL of EDTA plasma ambient in a plastic vial.

### Specimen Minimum Volume

0.75 mL Note: This volume does not allow for repeat testing.

### Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	NA

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)	7 days	
	Refrigerated	14 days	
	Frozen	14 days	

## Clinical & Interpretive

### Reference Values

Age	Range
<18 y	Not Established
Adult Males	10.7-32.0
Adult Females	6.8-29.3

Glycemic control goal for diabetic patients: >10

GlycoMark is intended for use with managing glycemic control in diabetic patients. A low result corresponds to high glucose peaks.

1, 5-AG blood levels can be affected by clinical conditions or medications.

## Performance

### PDF Report

No

### Day(s) Performed

Monday through Friday

### Report Available

5 to 12 days

### Performing Laboratory Location

Esoterix Endocrinology

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### CPT Code Information

84378

### LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FGLMA	GlycoMark	53835-5

---

Result ID	Test Result Name	Result LOINC® Value
Z0275	GlycoMark	53835-5