

Overview

**Method Name**  
Enzymatic, Colorimetric Assay

**NY State Available**  
Yes

Specimen

**Specimen Type**  
Varies

**Specimen Required**  
Submit only 1 of the following specimens:

- Serum**  
Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 1 mL of serum ambient in a plastic vial.
- Plasma**  
Draw blood in a lavender-top (EDTA) tube(s). Spin down and send 1 mL of EDTA plasma ambient in a plastic vial.

**Specimen Minimum Volume**  
0.75 mL Note: This volume does not allow for repeat testing.

Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)	7 days	
	Refrigerated	14 days	
	Frozen	14 days	

Clinical & Interpretive

Reference Values

Age	Range
<18 y	Not Established
Adult Males	10.7-32.0
Adult Females	6.8-29.3
Glycemic control goal for diabetic patients: >10	

GlycoMark is intended for use with managing glycemic control in diabetic patients. A low result corresponds to high glucose peaks.

1, 5-AG blood levels can be affected by clinical conditions or medications.

Performance

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

5 to 12 days

Performing Laboratory Location

Esoterix Endocrinology

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

84378

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FGLMA	GlycoMark	53835-5

Result ID	Test Result Name	Result LOINC® Value
Z0275	GlycoMark	53835-5