

Reporting Title: Arsenic/Creat w/Reflex, Random,U  
Performing Location: Rochester

Specimen Requirements:

Patient Preparation:

- 1. For the 48-hour period prior to start of collection, patient **should not** eat seafood.
- 2. High concentrations of gadolinium and iodine are known to interfere with most metal tests. If either gadolinium- or iodine-containing contrast media has been administered, **a specimen should not be collected for 96 hours.**

Supplies: Urine Tubes, 10 mL (T068)

Collection Container/Tube: Clean, plastic urine container with no metal cap or glued insert

Submission Container/Tube: Plastic, 10-mL urine tube or clean, plastic aliquot container with no metal cap or glued insert

Specimen Volume: 6 mL

Collection Instructions:

- 1. Collect urine a random urine specimen.
- 2. See [Metals Analysis Specimen Collection and Transport](#) for complete instructions.

Specimen Type	Temperature	Time	Special Container
Urine	Ambient	72 hours	
	Refrigerated (preferred)	28 days	
	Frozen	28 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
CRETR	Creatinine, Random, U	Numeric	mg/dL	2161-8
608900	Arsenic/Creatinine Ratio, U	Numeric	mcg/g Cr	13463-5
608901	Total Arsenic Concentration	Numeric	mcg/L	5586-3

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

Components:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
ASCU	Arsenic/Creatinine Ratio, U	1	82175	Yes	No
CRETR	Creatinine, Random, U	1	82570	Yes	No

CPT Code Information:

82175

82570

Reflex Tests:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
SPAS	Arsenic Speciation, Random, U	1	82175	No	Yes

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
SPAS	607692	Inorganic Arsenic (Toxic)	Alphanumeric	mcg/L	12481-8
SPAS	607693	Methylated Arsenic (Toxic)	Alphanumeric	mcg/L	53779-5
SPAS	607694	Organic Arsenic (Non-Toxic)	Alphanumeric	mcg/L	53778-7
SPAS	609405	Toxic Arsenic Concentration	Alphanumeric	mcg/L	96250-6
SPAS	607691	Total Arsenic Concentration	Numeric	mcg/L	5586-3

Reference Values:

ARSENIC/CREATININE:

0-17 years: Not established

> or =18 years: <24 mcg/g creatinine

CREATININE:

> or =18 years: 16-326 mg/dL

Reference values have not been established for patients who are younger than 18 years of age.