



Test Definition: CLAV

Cortisol, Left Adrenal Vein, Serum

Reporting Title: Cortisol, LAV

Performing Location: Mayo Clinic Laboratories - Rochester Superior Drive

Ordering Guidance:

For confirming the presence of synthetic steroids, order SGSS / Synthetic Glucocorticoid Screen, Serum.

Specimen Requirements:

Supplies: Sarstedt Aliquot Tube, 5 mL (T914)

Collection Container/Tube: Red top

Specimen Volume: 1.5 mL Serum

Submission Container/Tube: Plastic vial

Collection Instructions:

1. Morning (8 a.m.) and afternoon (4 p.m.) specimens are preferred.
2. Include time of collection
3. Centrifuge and aliquot serum into a plastic vial.

Additional Information: If multiple specimens are collected, send a separate order for each specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	28 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
6346	Cortisol, LAV	Numeric	mcg/dL	2143-6

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

82533

Reference Values:

No established reference values