

Test Definition: UBT

Helicobacter pylori Breath Test

Reporting Title: H. pylori C Urea Breath Test

Performing Location: Mayo Clinic Laboratories - Rochester Superior Drive

Ordering Guidance:

An alternative test for the diagnosis of active *Helicobacter pylori* infection in patients is the HPFRP / *Helicobacter pylori* with Clarithromycin Resistance Prediction, Molecular Detection, PCR, Feces, which requires a different collection.

Specimen Requirements:

Patient Preparation:

- 1. Do not administer this test if this list of instructions is not followed, as test results may be inaccurate:
- a. Do not administer this test if patient is allergic to citric acid. Note: Product contains aspartame.
- **b. For 2 weeks prior to testing, patients should not** take proton-pump inhibitors (eg, Prilosec, Prevacid, Aciphex, Protonix, and Nexium), histamine 2-receptor antagonists (H2RA), (eg Pepcid, Tagamet, Axid, or Zantac), or bismuth compounds (eg, Pepto-Bismol).
- c. For 4 weeks prior to testing, patients should not take antibiotics.
- 2. Carafate (sucralfate) does not interfere with this test.
- 3. Use of antacids does not affect the accuracy of this test.

Fasting: 1 hour, required

Supplies: H pylori Breath Kit - Meridian BreathID (T906; fees apply)

Collection Instructions:

- 1. Do not collect if patient is younger than 3 years.
- 2. Follow instructions included with kit.
- 3. Mixing the (13)C-Urea Tablet
- a. Dissolve the Citrica and the (13)C-enriched urea tablet in 150 to 200 mL (5.1 to 6.8 oz.) of tap water in the provided drinking cup.
- b. Close the lid firmly using both hands. Place fingers over lid and shake thoroughly for a few minutes, until the Citrica Powder and the urea tablet are completely dissolved.

Note: Tiny particles may remain visible after thorough mixing. However, if more substantial particulate matter is still present after five minutes of mixing, discard the solution and repeat the procedure with a new kit.

Forms:

If not ordering electronically, complete, print, and send 1 of the following forms with the specimen:

- -Infectious Disease Serology Test Request (T916)
- -Gastroenterology and Hepatology Test Request (T728)
- -Microbiology Test Request (T244)
- -General Request (T239)

Specimen Type	Temperature	Time	Special Container	
Breath	Ambient	14 days	BREATH TEST BAG	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
81590	H. pylori C Urea Breath Test	Alphanumeric		29891-9



Test Definition: UBT

Helicobacter pylori Breath Test

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

83013

Reference Values:

Negative

Reference values apply to all ages.