

Test Definition: FSTSC

Streptozyme Screen with Reflex to Titer

Reporting Title: Streptozyme Screen w/Reflex Titer

Performing Location: Focus Diagnostics,

Specimen Requirements:

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	
Serum	Refrigerated (preferred)	7 days	
	Frozen	30 days	
	Ambient	8 hours	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
Z2239	Streptozyme Screen	Alphanumeric		In Process

CPT Code: 1 × 86403

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
FSTTI	Streptozyme Titer	1	86406	No	No

Result Codes for Reflex Tests:

Test ID	et ID Result ID Reporting Name		Туре	Unit	LOINC®
FSTTI	Z2240	Streptozyme Titer	Alphanumeric		In Process



Test Definition: FSTSC

Streptozyme Screen with Reflex to Titer

Reference Values:

Negative

Test Performed by: Focus Diagnostics, Inc.

5785 Corporate Ave. Cypress, CA 90630-4750