

---

**Reporting Title:** Tobramycin, Trough, S  
**Performing Location:** Rochester**Specimen Requirements:**

Collection Container/Tube:

Preferred: Serum gel

Acceptable: Red top

Specimen Volume: 0.5 mL

Collection Instructions:

1. Draw blood no more than 30 minutes before next scheduled dose.
2. Serum gel tubes should be centrifuged within 2 hours of collection.
3. Red-top tubes should be centrifuged and aliquoted within 2 hours of collection.

Specimen Type	Temperature	Time
Serum	Refrigerated (preferred)	7 days
	Frozen	28 days
	Ambient	72 hours

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
TOBTA	Tobramycin, Trough, S	Numeric	mcg/mL	In Process

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:****CPT Code:** 1 × 80200**Reference Values:**

Therapeutic: &lt;2.0 mcg/mL

Toxic: &gt;2.0 mcg/mL