

TEST STATUS**NOTIFICATION DATE:** May 2, 2017**EFFECTIVE DATE:** June 2, 2017

CHLAMYDIA AND CHLAMYDOPHILA ANTIBODY PANEL 3
(IgG, IgA, IgM)
Test ID: CHPA

EXPLANATION: Effective June 2, 2017, test ID CHPA, referred to Quest Diagnostics Infectious Disease, Inc., will become obsolete due to the availability of Mayo test ID SCLAM, which includes IgM- and IgG-class antibodies to these bacteria. The clinical utility of evaluation for IgA-class antibodies is limited and not recommended by any national or international guidelines.

RECOMMENDED ALTERNATIVE TEST:

CHLAMYDIA SEROLOGY, SERUM
Test ID: SCLAM

METHODOLOGY: Micro-Immunofluorescent Antibody (MIF) Assay

REFERENCE VALUES:

<i>Chlamydomphila pneumoniae</i>	<i>Chlamydomphila psittaci</i>	<i>Chlamydia trachomatis</i>
IgG: <1:64	IgG: <1:64	IgG: <1:64
IgM: <1:10	IgM: <1:10	IgM: <1:10

SPECIMEN REQUIREMENTS:**Container/Tube:****Preferred:** Serum gel**Acceptable:** Red top**Specimen Volume:** 0.2 mL**Minimum Volume:** 0.15 mL**SPECIMEN STABILITY INFORMATION:**

Specimen Type	Temperature	Time
Serum	Refrigerated (preferred)	30 days
	Frozen	30 days

CPT CODE: 86631 x 3-IgG
86632 x 3-IgM

DAY(S) SET UP: Monday through Saturday **ANALYTIC TIME:** Same day/1 day

QUESTIONS: Contact Mary Erath, MML Referrals Supervisor
Telephone: 800-533-1710