

## Motilin, Plasma or Serum

**Test ID:** FMOT

**Explanation:**

Effective August 15, 2019 Test ID FMOT, referred to Inter Science Institute, will become obsolete.

**Recommended Alternative Test:**

## Miscellaneous Inter Science Institute

**Test ID:** ZW71

**Referral Lab Code:** MML1447

**Test Name:** Motilin (Plasma/Serum)

**Methodology:**

Direct Radioimmunoassay (RIA)

**Reference Values:**

Interpretative report

**Specimen Requirements:**

**Patient preparation:** Patient should be fasting 10-12 hours prior to collection of specimen. Patient should not be on any antacid medication or medications affecting intestinal motility, if possible, for at least 48 hours prior to collection.

**Container/Tube:** Red-top, serum gel tube or Lavender-top EDTA tube

**Preferred:** Serum or Plasma

**Acceptable:** NA

**Specimen Volume:** 3 mL

**Collection Instructions:** Serum: Draw blood in a plain red-top or serum gel tube(s). Separate immediately and send 3 mL of serum frozen in a plastic vial.  
Plasma: Draw blood in a lavender-top (EDTA) tube(s). Separate immediately and send 3 mL of EDTA plasma frozen in a plastic vial.

**Minimum Volume:** 1 mL

**Specimen Stability Information:**

Specimen Type	Temperature	Time
Varies	Frozen	90 days
	Refrigerated	7 days

**CPT Code:**

83519

**Miscellaneous Referral Test Fee:**

\$235.00

**Day(s) Setup:** Monday-Friday    **Analytic Time:** <28 days**Questions:**

Contact Mary Erath, MCL Referrals Supervisor at 800-533-1710.