

Notification Date: July 22, 2019 Effective Date: August 15, 2019

## Ganciclovir, Serum

Test ID: GANC

## **Explanation:**

Due to low test utilization, test code GANC will become obsolete on the effective date.

#### **Recommended Alternative Test:**

# Ganciclovir, Plasma

**Test ID: FGANP** 

## Methodology:

High Performance Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS)

### **Reference Values:**

0.2 mcg/mL

The reported value may vary based on dilution or weight of the specimen

#### Specimen Requirements:

Container/Tube: EDTA (lavender top) or Pink top

Specimen Volume: 3 mL

Draw blood in EDTA (lavender top) or Pink top tube(s) (Plasma gel tube is not acceptable.) Spin down and

send 1 mL of plasma refrigerate in preservative-free plastic vial.

## **Specimen Stability Information:**

| Specimen Type | Temperature              | Time     |
|---------------|--------------------------|----------|
| Plasma EDTA   | Refrigerated (preferred) | 30 days  |
|               | Ambient                  | 30 days  |
|               | Frozen                   | 120 days |

#### **CPT Code:**

80375

**Day(s) Setup:** Monday – Sunday **Analytic Time:** 8 days

Questions: Contact Kim Terrio, Laboratory Technologist Resource Coordinator at 800-533-1710.