

Reporting Title: BRCANext-Expanded**Performing Location:** Ambry Genetics**Specimen Requirements:**

Container/Tube:

Preferred: EDTA (purple top)

Acceptable: Yellow top (citric acetate) or grey top (potassium oxalate/sodium fluoride) tube

Specimen volume: 10 mL

Collection Instructions: Draw blood in EDTA (purple top) tube(s) and send 10 mL of whole blood refrigerate.

Complete Ambry Cancer test Comprehensive requisition form.

NOTE: For transfusion patients, wait at least 2 weeks after a packed cell or platelet transfusion and at least 4 weeks after a whole blood transfusion prior to blood draw.

Specimen Type	Temperature	Time	Special Container
Whole Blood EDTA	Refrigerated (preferred)	7 days	
	Ambient	72 hours	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
FBRAC	BRCANext-Expanded	Alphanumeric		Not Provided

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

Referral

CPT Code: 1 × 81432