

Glucose, Random, Urine

Test ID: RGLUR

Explanation:

On the effective date this test will be replaced by GLUR1/ Glucose, Random, Urine.

Recommended Alternative Test:

Glucose, Random, Urine

Test ID: GLUR1

Methodology:

Reflectance Spectrophotometry

Reference Values:

Negative

Specimen Requirements:

Container/Tube:	Plastic, 5-mL tube
Specimen Volume:	5 mL
Collection Instructions:	1. Collect a random urine specimen. 2. No preservative
Minimum Volume:	2 mL

Specimen Stability Information:

Specimen Type	Temperature	Time
Urine	Refrigerated (preferred)	72 hours
	Ambient	2 hours

CPT Code:

82945

Day(s) Setup: Monday through Sunday; Continuously **Analytic Time:** Same day/1 day

Questions

Contact Nancy Benson, Laboratory Technologist Resource Coordinator at 800-533-1710.