
Diuretic Screen, Urine

Test ID: FDIRU**Explanation:**

Effective August 25, 2025, Test ID FDIRU, referred to Medtox Laboratory, will be obsolete and replaced with new test FDISU due to file definition change.

Recommended Alternative Test:

Diuretic Profile, Urine

Test ID: FDISU**Methodology:**

High Performance Liquid Chromatography with Ultraviolet Detection; (HPLC-UV)

Specimen Requirements:**Specimen Type:** Urine**Container/Tube:** Plastic preservative-free urine container**Specimen Volume:** 10 mL**Collection Instructions:**

1. Collect 10 mL random urine without preservative.
2. Send specimen refrigerated

Specimen Stability Information:

Specimen Type	Temperature	Time
Urine	Refrigerated (preferred)	14 days
	Ambient	72 hours
	Frozen	14 days

CPT Code:

80377

Day(s) Setup: Thursday**Report Available:** 7 to 9 days**Questions**

Contact MCS Referrals Supervisors Amy Bluhm or Heather Hilgart at 800-533-1710.