# MAYO CLINIC LABORATORIES

### **TEST STATUS - TEST DOWN**

Notification Date: October 8, 2025 Effective Date: Immediately

# Temazepam and Oxazepam, Serum

Test ID: TEMZS

**Explanation:** This test will made non-orderable effective immediately due to unexpected analytic issues. An update will be provided when testing resumes.

#### **Recommended Alternative Test:**

## Temazepam, serum or plasma

Test ID: ZW80

Referral Lab Code: 155

### Methodology:

Liquid chromatography/tandem mass spectrometry (LC/MS-MS)

## Specimen Requirements:

Specimen Type: Serum

Collection Container/Tube: Red Top

Submission Container/Tube: Preservative-Free plastic vial

Specimen Volume: 2.0 mL Minimum Volume: 0.5 mL Collection Instructions:

1. Separate and transfer serum into a plastic screw capped vial.

2. Ship Refrigerated

### Specimen Stability Information:

Specimen Type	Temperature	Time
Serum Red	Refrigerate (preferred)	7 days
	Ambient	72 hours
	Frozen	180 days

**CPT Code:** 80299

**Miscellaneous Referral Test Fee:** \$69.57

**Day(s) Performed:** Monday – Sunday Report Available: 5 – 7 days

### Questions

Contact Brandon DeBoom, Laboratory Resource Coordinator at 800-533-1710.