

Amphetamine-Type Stimulants Confirmation, Chain of Custody, Meconium

Test ID: AMPMX

Explanation: Due to method enhancements, cutoff concentrations for analytes reported within this test will be lowered from 50 ng/g to 20 ng/g.

Current Reference Value	New Reference Value
Negative Positives are reported with a quantitative LC-MS/MS result. Cutoff concentrations: Amphetamine by LC-MS/MS: 50 ng/g Methamphetamine by LC-MS/MS: 50 ng/g 3,4-Methylenedioxyamphetamine by LC-MS/MS: 50 ng/g 3,4-Methylenedioxyethylamphetamine by LC-MS/MS: 50 ng/g 3,4-Methylenedioxymethamphetamine by LC-MS/MS: 50 ng/g	Negative Positives are reported with a quantitative liquid chromatography-tandem mass spectrometry (LC-MS/MS) result. Cutoff concentrations for LC-MS/MS testing: Amphetamine: 20 ng/g Methamphetamine: 20 ng/g 3,4-Methylenedioxyamphetamine: 20 ng/g 3,4-Methylenedioxyethylamphetamine: 20 ng/g 3,4-Methylenedioxymethamphetamine: 20 ng/g

Questions

Contact Brandon DeBoom, Laboratory Technologist Resource Coordinator at 800-533-1710.