

## Immunoglobulin G (IgG), Serum

**Test ID:** IGG

**Explanation:** Ambient temperature specimen stability will be updated to 14 days on the effective date.

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	28 days
	Ambient	28 days
	Frozen	28 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	28 days
	Ambient	14 days
	Frozen	28 days

### Questions

Contact Michelle Rath, Laboratory Technologist Resource Coordinator at 800-533-1710.