

Inhibin A, Tumor Marker, S

Test ID: INHA

Explanation: Updating specimen type and stability as indicated below.

Current Specimen Required	
Collection Container/Tube:	
Preferred:	Red Top
Acceptable:	Serum Gel

New Specimen Required	
Collection Container/Tube:	
Preferred:	Serum Gel
Acceptable:	Red Top

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	7 days
	Frozen	90 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	14 days
	Frozen	90 days

Questions

Contact Bethany Feind, Laboratory Technologist Resource Coordinator at 800-533-1710.