

Urine Drugs of Abuse and Compliance Monitoring Testing

Explanation: As a result of pre-analytic process improvements and standardization efforts, the specimen volume requirements will be changing for the list of tests below. The following specimen submission containers are available in our supply catalog if needed:

T313 – Urine Container, 60 mL

T068 – Urine Tubes, 10 mL

T465 – Aliquot Tube, 5 mL

Mayo Test ID	Reporting Name	Current Requested Volume	New Requested Volume	Current Minimum Volume	New Minimum Volume
ADLTX	Adulterants Survey, CoC, U	20.0 mL	1.5 mL	0.4 mL	1.0 mL
ADULT	Adulterants Survey, U	20.0 mL	1.5 mL	0.4 mL	1.0 mL
AMPHU	Amphetamines Confirmation, U	5.0 mL	1.0 mL	2.5 mL	0.5 mL
AMPHX	Amphetamines Confirmation, CoC, U	20.0 mL	5.0 mL	2.5 mL	1.0 mL
BARBX	Barbiturates Confirmation, CoC, U	20.0 mL	10.0 mL	5.0 mL	5.0 mL
BNZU	Benzodiazepines Confirmation, U	10.0 mL	1.0 mL	5.0 mL	0.5 mL
BNZX	Benzodiazepines Conf, CoC, U	20.0 mL	5.0 mL	5.0 mL	1.0 mL
BUPM	Buprenorphine and Metabolite, U	5.0 mL	1.0 mL	1.0 mL	0.5 mL
BUPMX	Buprenorphine w/metabolite, CoC, U	20.0 mL	5.0 mL	1.0 mL	1.0 mL
COKEX	Cocaine and metabolite Conf, CoC, U	20.0 mL	10.0 mL	5.0 mL	5.0 mL
ETGC	Ethyl Glucuronide Confirmation, Random, Urine	5.0 mL	1.0 mL	2.5 mL	0.5 mL
ETGR	Ethyl Glucuronide Scrn w/Reflex, U	20.0 mL	5.0 mL	2.5 mL	2.5 mL
ETGS	Ethyl Glucuronide Screen, U	20.0 mL	2.0 mL	2.5 mL	0.5 mL
ETGX	Ethyl Glucuronide Conf, CoC, U	20.0 mL	5.0 mL	2.5 mL	1.0 mL
FENR	Fentanyl Screen w/Reflex, U	20.0 mL	5.0 mL	2.5 mL	2.5 mL
FENS	Fentanyl Screen, U	20.0 mL	2.0 mL	2.5 mL	0.5 mL
FENTU	Fentanyl w/metabolite Conf, U	3.0 mL	10.0 mL	2.1 mL	2.1 mL
FENTX	Fentanyl w/metabolite Conf, CoC, U	20.0 mL	10.0 mL	2.1 mL	2.1 mL
HYDCU	Hydrocodone with Metabolite Confirmation, Random, Urine	20.0 mL	3.0 mL	2.5 mL	1.0 mL
HYDMU	Hydromorphone Confirmation, Random, Urine	20.0 mL	3.0 mL	2.5 mL	1.0 mL

MTDNX	Methadone Confirmation, CoC, U	20.0 mL	10.0 mL	2.5 mL	2.5 mL
OPATU	Opiate Confirmation, U	3.0 mL	3.0 mL	2.5 mL	1.0 mL
OPATX	Opiate Confirmation, CoC, U	20.0 mL	5.0 mL	2.5 mL	2.0 mL
OXYCU	Oxycodone with Metabolite Confirmation, Random, Urine	20.0 mL	3.0 mL	2.5 mL	1.0 mL
OXYCX	Oxycodone w/metabolite Conf, CoC, U	20.0 mL	5.0 mL	2.5 mL	2.0 mL
OXYMU	Oxymorphone Confirmation, U	20.0 mL	3.0 mL	2.5 mL	1.0 mL
OXYSU	Oxycodone Screen, U	20.0 mL	2.0 mL	2.5 mL	0.5 mL
TAPEN	Tapentadol and Metabolite, U	2.0 mL	1.0 mL	0.1 mL	0.5 mL
TBSU	Targeted Benzodiazepine Screen, U	3.0 mL	1.0 mL	1.0 mL	0.5 mL
TOSU	Targeted Opioid Screen, U	3.0 mL	1.0 mL	1.0 mL	0.5 mL
TRAM	Tramadol and Metabolite, U	2.0 mL	1.0 mL	0.1 mL	0.5 mL
TSPU	Targeted Stimulant Screen, U	3.0 mL	1.0 mL	1.0 mL	0.5 mL

Questions

Contact Brandon DeBoom, Laboratory Technologist Resource Coordinator at 800-533-1710