

Aspergillus Antigen, Bronchoalveolar Lavage

Test ID: ASPBA

Explanation: On effective date, preferred specimen stability will be updated to refrigerate.

Current Specimen Stability		
Specimen	Temperature	Time
Lavage	Frozen (preferred)	14 days
	Refrigerated	14 days

New Specimen Stability		
Specimen	Temperature	Time
Lavage	Refrigerated (preferred)	14 days
	Frozen	14 days

Questions

Contact Dunisha Messmer Laboratory Technologist Resource Coordinator at 800-533-1710.