

## Hematologic Disorders, Fluorescence In Situ Hybridization (FISH) Hold, Varies

**Test ID:** HOLDF

**Explanation:**

On the effective date the preferred specimen container/tube for all accepted specimen types will be updated as shown below.

Current Specimen Container/Tube Required
Green top (sodium heparin)

New Specimen Container/Tube Required	
Preferred	Yellow top (ACD)
Acceptable	Green top (sodium heparin)
Acceptable	Lavender top (EDTA)

**Questions**

Contact Joshua Couchene Laboratory Technologist Resource Coordinator at 800-533-1710.