## **TEST STATUS**



Notification Date: December 6, 2021 Effective Date: Immediately

# Immunoglobulin D (IgD), Serum

Test ID: IGD

## **Explanation:**

Due to a reagent shortage, IGD: Immunoglobulin D (IgD), Serum will become nonorderable effective immediately. A notification will be sent when testing resumes.

#### **Recommended Alternative Test:**

# Immunoglobulin D, Serum

Test ID: FIMMD

### Methodology:

Quantitative Immunoturbidimetry

#### **Reference Values:**

Less than or equal to 15.3 mg/dL

## **Specimen Requirements:**

Preferred: SST

Acceptable: Red-top tube

Specimen Volume: 1 mL

Minimum Volume: 0.5 mL

Collection Instructions: Draw blood in a serum gel tube(s), red-top tube(s) is acceptable. Separate serum

from cells within 2 hours of collection. Transfer 1 mL of serum into a plastic vial.

Ship refrigerated.

## **Specimen Stability Information:**

Specimen Type	Temperature	Time		
Serum SST	Refrigerated (preferred)	14 days		
	Frozen	180 day		

_	_	 _				
"	P٦	 ~	o	М	Δ	•
v		•	v	ч	◡	

82784

## **Miscellaneous Referral Test Fee:**

\$87.70

Day(s) Performed: Monday, Wednesday, Friday Report Available: 1 to 12 days

# **Questions**

Contact Amy Ennis, Laboratory Technologist Resource Coordinator at 800-533-1710.