

## Alpha-Fetoprotein (AFP), Single Marker Screen, Maternal, Serum

Test ID: MAFP1

**Explanation:** Updating specimen type, stability as indicated below.

Current Specimen Required	
Preferred	Red top
Acceptable	Serum Gel

New Specimen Required	
Preferred	Serum gel
Acceptable	Red top

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	7 days
	Frozen	90 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	14 days
	Ambient	7 days
	Frozen	90 days

### Questions

Contact Bethany Feind, Laboratory Technologist Resource Coordinator at 800-533-1710.