

Notification Date: February 16, 2022

Effective Date: Immediate

Aripiprazole (Abilify)

Test ID: FARI

Explanation: The following changes for test ID FARI, sent to Medtox, are effective immediately.

Current Specimen Stability			
Specimen	Temperature	Time	
Varies	Refrigerated (preferred)	7 days	
	Ambient	72 hours	
	Frozen	180 days	

New Specimen Stability			
Specimen	Temperature	Time	
Varies	Refrigerated (preferred)	7 days	
	Ambient	72 hours	
	Frozen	25 days	

Current CPT Code	
80342	

New CPT Code	
80299	

Questions:

Contact MCL Referrals Supervisor Amy Bluhm at 800-533-1710.