

Greiner® Z Trace Element no-additive (Aluminum Only), 6 mL

T#: T713

Explanation:

Due to a national backorder, we are temporarily out of stock item T713, the Greiner Z Trace Element no-additive tube, used for Aluminum testing. There are no alternate tubes available.

Tests:

Aluminum, Serum

Reference Values:

<7 ng/mL

<60 ng/mL (dialysis patients)

Note:

The Royal Blue (BD) top tubes should NOT be used for this test. Specimens are contaminated by the stopper.

Questions

Contact Rebekah Knauer, Laboratory Technologist Resource Coordinator at 800-533-1710.