

NEW REFERRAL TEST

Notification Date: March 30, 2022 Effective Date: April 15, 2022

Hypersensitivity Pneumonitis Panel

Explanation: The following new orderable, referred to Eurofins Viracor, will be available on 04/15/2022.

Test ID: FHSPP

Methods:

Immunoassay FEIA

Reference Values:

Alternaria tenuis/alternata IgG <12.0 mcg/mL Aspergillus fumigatus IgG <46.0 mcg/mL Aureobasidium pullulans IgG <18.0 mcg/mL Laceyella sacchari IgG <13.0 mcg/mL Micropolyspora faeni IgG <5.0 mcg/mL Penicillium Chrysogenum/notatum IgG <22.0 mcg/mL Phoma betae IgG <8.0 mcg/mL Trichoderma viride IgG <10.0 mcg/mL

0.6 mL

Specimen Requirements:

Container/Tube:

Specimen Volume:

Preferred: Red top

Acceptable: Serum gel

Collection Instructions: Draw blood in a plain, red-top tube(s), serum gel tube(s) is acceptable. Spin down

and send 0.6 mL of serum refrigerated in a plastic vial.

Minimum Volume: 0.5 mL

Specimen Stability Information:

Specimen Type	Temperature	Time
Serum	Refrigerated (preferred)	28 days
	Ambient	7 days
	Frozen	28 days

Fee: CPT Code: \$153.40 86001 x8

Day(s) Setup: Monday through Friday Report Available: 2 to 7 days

Questions:

Contact MCL Referrals Supervisor Amy Bluhm at 800-533-1710.