

Nicotine and Metabolites, Random, Urine

Test ID: NICOU

Explanation: On the effective date, the preferred specimen temperature for NICOU will be changed to refrigerated.

Current Specimen Stability		
Specimen	Temperature	Time
Urine	Ambient (preferred)	28 days
	Refrigerated	28 days
	Frozen	365 days

New Specimen Stability		
Specimen	Temperature	Time
Urine	Refrigerated (preferred)	28 days
	Ambient	28 days
	Frozen	365 days

Questions

Contact Joshua Yang, Laboratory Technologist Resource Coordinator at 800-533-1710.