

## Hypersensitivity Pneumonitis Avian Panel

**Test ID:** FAVI

**Explanation:**

Effective immediately, test ID FAVI, referred to Medical College of WI-Allergy-Immunology Diagnostic Lab, will reflect the following changes:

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	365 days
	Ambient	24 hours
	Frozen	365 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (Preferred)	30 days
	Frozen	365 days

Current Days Performed
Monday through Friday

New Days Performed
Thursday and Friday

Current Analytic Time
8 Days

New Analytic Time
10-14 Days

**Questions:**

Contact MCL Referrals Supervisor Amy Bluhm at 800-533-1710.