

## ST2, Serum

Test ID: ST2S

**Explanation:** This test will be made non-orderable effective immediately due to reagent supply issues. There are no referral options available through Mayo Clinic Laboratories. Though an equivalent test is not available, GAL3 is listed below to consider as an alternative. An update will be provided when testing resumes.

### Suggested Alternative Test:

## Galectin-3, Serum

Test ID: GAL3

### Methodology:

Enzyme-Linked Immunosorbent Assay (ELISA)

### Reference Values:

<24 months: not established

2-17 years: < or =25.0 ng/mL

> or =18 years: < or =22.1 ng/mL

### Specimen Requirements:

|                            |   |
|----------------------------|---|
| Collection Container/Tube: | Red top   |
| Submission Container/Tube: | Plastic vial                                    |
| Specimen Volume:           | 1 mL  |
| Collection Instructions:   | Centrifuge and aliquot serum into plastic vial. |
| Minimum Volume:            | 0.2 mL  |

**Specimen Stability Information:**

| Specimen Type | Temperature        | Time     |
|---------------|--------------------|----------|
| Serum Red     | Frozen (preferred) | 365 days |
|               | Refrigerated       | 24 hours |

**CPT Code:** 82777**Day(s) Performed:** Monday**Report Available:** 1 – 8 days**Questions**

Contact Nancy Benson, Laboratory Resource Coordinator at 800-533-1710.