MAYO CLINIC LABORATORIES

TEST STATUS - DOWN

Notification Date: February 15, 2023 Effective Date: Immediately

F12 Gene, Next-Generation Sequencing, Varies

Test ID: F12NG

Explanation:

Due to reagent supply issue, F12NG/F12 Gene, Next-Generation Sequencing, Varies will become nonorderable immediately. A notification will be sent once testing resumes.

Recommended Alternative Test:

PGxome Custom Panel

Test ID: ZW194

Referral Lab Code: 6000

NYS Approved: No

Methodology:

Copy Number Variant (CNV) Analysis via NGS data

Specimen Requirements:

Specimen Type: Whole blood

Container/Tube:

Preferred: Lavender top (EDTA)

Acceptable: ACD

Specimen Volume: 3-5mL

Minimum Volume: 3ml, (Small Infants: 1ml)

Note:

Complete and submit with sample: Prevention Genetics Lab Request Form

Specimen Stability Information:

Specimen Type	Temperature	Time
Whole Blood	Refrigerated	7 days

CPT Code:

Call for CPT code(s)

Miscellaneous Referral Test Fee:

Call for pricing

Day(s) Performed: Varies Report Available: 25 days

Questions

Contact Connie Penz, Laboratory Resource Coordinator at 800-533-1710.