

Thyrotropin Receptor Antibody, Serum

Test ID: THYRO

Explanation: Updating stability, specimen type etc as indicated below.

Current Specimen Required	
Preferred	Red top
Acceptable	Serum Gel

New Specimen Required	
Preferred	Serum Gel
Acceptable	Red top

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	7 days
	Frozen	30 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	7 days
	Frozen	90 days

Questions

Contact Bethany Feind, Laboratory Resource Coordinator at 800-533-1710.