

Prenatal Hepatitis Evaluation, Serum

Test ID: PHSP

Explanation: Due to a method change, the following changes will be implemented on the effective date.

Current Method
Chemiluminescence Immunoassay (CIA)

New Method
Electrochemiluminescence Immunoassay (ECLIA)

Current Specimen Stability		
Specimen	Temperature	Time
Serum SST	Frozen (preferred)	28 days
	Ambient	24 hours
	Refrigerated	7 days

New Specimen Stability		
Specimen	Temperature	Time
Serum SST	Frozen (preferred)	84 days
	Refrigerated	6 days

Questions

Contact Dunisha Messmer, Laboratory Resource Coordinator at 800-533-1710.